
## Bedienungsanleitung VET 112 Digitales Veterinär-Thermometer Kat. Nr. 15.2019

# TFA®

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### 1. Before you use this product

**Please make sure you read the instruction manual carefully.**

- Following and respecting the instructions in your manual will prevent damage to your instrument and loss of your statutory rights arising from defects due to incorrect use.
- We shall not be liable for any damage occurring as a result of non following of these instructions. Likewise, we take no responsibility for any incorrect readings or for any consequences resulting from them.
- Please take particular note of the safety advice!
- Please keep this instruction manual safe for future reference.

### 2. For your safety

- This product is intended exclusively for rectal measurement of the body temperature of large animals. It should only be used as described within these instructions.
- Check the thermometer for any signs of damage or wear before use. Use the unit only when it is perfectly fine.
- Unauthorized repairs, alterations or changes to the product are prohibited.
- Operation in strong electromagnetic fields, such as next to a mobile phone, can lead to malfunction. Keep a minimum distance of 3.3 m from such devices when using this device.
- Always discuss the determined values with your veterinarian.

**Caution!**  
**Risk of injury:**

- The instrument contains small parts that could be swallowed. Keep this device and the battery out of reach of children and animals.
- Batteries contain harmful acids and may be hazardous if swallowed. If a battery is swallowed, this can lead to serious internal burns and death within two hours. If you suspect a battery could have been swallowed or otherwise caught in the body, seek medical help immediately.
- Batteries must not be thrown into a fire, short-circuited, taken apart or recharged. Risk of explosion!
- Low batteries should be changed as soon as possible to prevent damage caused by leaking.
- Wear chemical-resistant protective gloves and safety glasses when handling leaking batteries.

**Important information on product safety!**

- Do not expose the device to extreme temperatures, vibrations or shocks. Do not bend the instrument.

### 3. Elements and buttons

- A:** Battery compartment  
**C:** Display  
**C 2:** Battery symbol  
**D:** Temperature probe
- B:** ON/OFF button  
**C 1:** Temperature display  
**C 3:** Symbol measurement memory

### 4. Use


- Disinfect the tip before each use!
- Press the ON/OFF button.
- When the unit is switched on, a signal tone sounds and all segments are briefly displayed. The instrument starts a self-test.
- The latest measured temperature (M) will appear shortly on the display, then - - - °C.
- If the room temperature is higher than 32 °C, the room temperature would be displayed instead of - - - °C. °C starts flashing. The device is now ready for use.
- Place the thermometer at the site of measurement.
- The actually measured maximum temperature is permanently shown on the display and °C will be flashing.
- If the maximum temperature remains stable, °C is no longer flashing and an alarm tone will sound.
- Remove the thermometer. Press the ON/OFF button for 3 seconds to switch off the device. The device will turn off automatically after approx. 9 minutes.

### 5. Care and maintenance

- Clean the instrument before and after each use with a soft damp cloth. For disinfection please use water-diluted isopropyl alcohol according to the manufacturer's instructions. Any other cleansing substances or cleansing method can cause dysfunction or damage to the instrument. Do not immerse in water or sterilize by boiling method, gas or steam autoclaves.

**Instruction manual**  
**VET 112**  
**Digital veterinary thermometer**  
**Art.-No. 15.2019**

**TFA**

- Remove the battery if you do not use the device for a long period of time.
- Keep the device in a dry place.

### 5.1 Battery replacement

- Change the battery when the display starts to fade or the battery symbol "▼" appears on the display.
- Open the battery compartment by turning the cover to the left using a coin (until you reach the upper point).
- Remove the battery and insert a new battery LR 44 (+ pole above). Make sure the polarity is correct.
- Close the battery cover again.

### 6. Troubleshooting

Problem	Solution
No display	→ Switch on the device → Ensure that the battery polarity is correct (+pole above) → Change the battery
Incorrect indication	→ Change the battery → Check the correct position of the thermometer
Display "Err"	→ Remove the battery Wait one minute and then reinsert the battery

### 7. Waste disposal

Never dispose of empty batteries and rechargeable batteries with ordinary household waste. As a consumer, you are required by law to take them to your retail store or to an appropriate collection site depending on national or local regulations in order to protect the environment.

The symbols for the contained heavy metals are: Cd = cadmium, Hg = mercury, Pb = lead

This product is labelled in accordance with the EU Waste Electrical and Electronic Equipment Directive (WEEE). Please do not dispose of this product in ordinary household waste.

As a consumer, you are required to take end-of-life devices to a designated collection point for the disposal of electrical and electronic equipment, in order to ensure environmentally-compatible disposal.

### 9. Specifications

Type	Maximum thermometer
Measuring range	32 °C ... 42.9 °C
--- °C	Surrounding temperature is lower than 32 °C
H °C	Surrounding temperature is higher than 42.9 °C
Accuracy	±0,1 °C at 35,5 °C...42 °C ±0,2 °C otherwise
Resolution	0.1 °C
Measurement time (rectal)	Approx. 20-30 seconds
Storage conditions	Temperature -25 °C...+60 °C humidity 15...95 % rH
Ambient conditions by use	Temperature 10 °C...+40 °C humidity 15...95 % rH
IP protection class	IP67
Battery	1 x LR44 1,5 V (included)
Expected service life	5 years or 10000 measurements
Housing dimension	100 x 20 x 195 mm
Weight	45 g

TFA Dostmann GmbH & Co.KG, Zum Ottersberg 12, D-97877 Wertheim, Germany  
No part of this manual may be reproduced without written consent of TFA Dostmann. The technical data are correct at the time of going to print and may change without prior notice.  
The latest technical data and information about this product can be found in our homepage by simply entering the product number in the search box.

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Fabricante: Microlife Corporation, 9F, 431 RuiGuang Road, NeiHu Taipei 11492, Taiwan, R.O.C.
Fabricado para: TFA Dostmann GmbH & Co. KG, Zum Ottersberg 12, 97877 Wertheim, Germany
TFA Dostmann GmbH & Co.KG, Zum Ottersberg 12, D-97877 Wertheim, Alemania