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Kat. Nr. 50.1003

TFA

SWING - Personenwaage

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TFA

SWING - Bathroom scale

Thank you for choosing this instrument from the firm TFA.

1. Before you use this product

- Please make sure you read the instruction manual carefully.
- Following and respecting the instructions in your manual will prevent damage to your instrument and loss of your statutory rights arising from defects due to incorrect use.
- We shall not be liable for any damage occurring as a result of non-following of these instructions. Likewise, we take no responsibility for any incorrect readings or for any consequences resulting from them.
- Please take particular note of the safety advice!
- Please keep this instruction manual safe for future reference.

2. Delivery contents

- Bathroom scale
- Batteries 3 x 1.5 V AA
- Instruction manual

3. Field of operation and all the benefits of your new instrument at a glance

- Analogue and digital display in tachometer style
- Automatic power-on function – no tapping required
- Automatic power-off function
- Modern flat design
- Tempered safety-glass surface
- Weighing reminder function with acoustic signal
- Backlight

4. For your safety

- This product is exclusively intended for the field of application described above. It should only be used as described within these instructions.
- Unauthorised repairs, modifications or changes to the product are prohibited.
- Before starting a treatment or diet discuss with your doctor first.

Caution! Risk of injury:

- Do not stand on the scale with wet feet and especially not when the surface of the scale is wet. There is a dangerous risk of slipping.
- Do not step on the edge of the scale: Danger of tilting!
- Keep the product and the batteries out of the reach of children.
- Batteries contain harmful acids and may be hazardous if swallowed. If a battery is swallowed, this can lead to serious internal burns and death within two hours. If you suspect a battery could have been swallowed or otherwise caught in the body, seek medical help immediately.
- Batteries must not be thrown into a fire, short-circuited, taken apart or recharged. Risk of explosion!
- Low batteries should be changed as soon as possible to prevent damage caused by leaking. Never use a combination of old and new batteries together, nor batteries of different types.
- Wear chemical-resistant protective gloves and safety glasses when handling leaking batteries.

Important information on product safety!

- Do not expose the device to extreme temperatures, vibrations or shocks.
- Do not place any objects on the unit.
- Protect it from humidity.

5. Set up and operation

- Open the battery compartment at the back and insert the three batteries 1.5 V AA with the polarity as indicated.
- Your scale is now ready to use.
- 88:88.8 appears on the display.
- Press the KG/LB button at the back and select your desired weight unit.
- Place the scale on a flat solid surface.
- In order to switch on the scale, just stand on the platform.
- The backlight turns on.
- After a few seconds your measured weight will be displayed in kg (default) or in the weighting unit you have set.
- The measuring process is finished as soon as your weight flashes on the display.
- The weight remains on the display after leaving the scale.
- The scale will turn off automatically after approximately 5 seconds.

5.1 Setting of the weight unit

- At the back you will find the KG/LB button to convert the weight unit. You can select between "kg" for Kilogram, "lb" for British and American pound (1 pound = 453,59237 Gram).

Instruction manual

Kat. Nr. 50.1003

TFA

5.2 Setting of the weighing reminder function

- At the back you will find the 24H/ALARM button to set the weighing reminder function.
- Press the button to activated the function
- On the display appears the alarm symbol, the backlight turns on shortly and a beep sounds.
- After 24 hours a beep (20 seconds) sounds.
- This ensures that for example during a diet, the weight is always checked at the same time and the result is not distorted by the normal weight fluctuations over the course of a day.
- Press the button again to deactivate the function.
- The alarm symbol will disappear.

5.3 Use of the scale

- If possible, always weigh yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results.
- When you start using the device and after replacing the batteries, the device has to be adjusted.
- To end the adjustment, repeat the procedure three times.
- After the adjustment on the display of the balance will appear after each new use the correct weight.

6. Care and maintenance

- Clean your instrument with a soft damp cloth. Do not use solvents or scouring agents.
- Remove the batteries if you do not use it for a long period of time.
- Keep your scale in a dry place.

7. Troubleshooting

Problems	Solution
No indication	→ Stand on the platform Ensure that the batteries' polarity is correct → Change the batteries
Indication "---"	→ Overload (max. 180 kg)
Indication "Lo"	→ Change the batteries
Incorrect display	→ Change the batteries

If your device fails to work despite these measures contact the supplier from whom you purchased it.

8. Waste disposal

This product has been manufactured using high-grade materials and components which can be recycled and reused.

Never dispose of empty and rechargeable batteries in household waste. As a consumer, you are legally required to take them to your retail store or to an appropriate collection site depending on national or local regulations in order to protect the environment.

The symbols for the heavy metals contained are: Cd=cadmium, Hg=mercury, Pb=lead

This instrument is labelled in accordance with the EU Waste Electrical and Electronic Equipment Directive (WEEE).

Please do not dispose of this instrument in household waste. The user is obligated to take end-of-life devices to a designated collection point for the disposal of electrical and electronic equipment, in order to ensure environmentally-compatible disposal.

9. Specifications

Measuring range:	5 – 180 kg / 11 - 396 lb
Resolution:	0.1 kg / 0.2 lb
Power consumption:	Batteries 3 x 1.5 V AA (included)
Dimensions:	315 x 32 x 315 mm
Weight:	1729 g (instrument only)

TFA Dostmann GmbH & Co. KG,
Zum Ottersberg 12, D-97877 Wertheim, Germany
No part of this manual may be reproduced without written consent of TFA Dostmann. The technical data are correct at the time of going to print and may change without prior notice.
The latest technical data and information about this product can be found in our homepage by simply entering the product number in the search box.
www.tfa-dostmann.de

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