Bedienungsanleitung
Kat. Nr. 50.1015.01

STEP PLUS – Personenwaage

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www.tfa-dostmann.de/en/service/downloads/instruction-manuals

STEP PLUS – Bathroom scale

Thank you for choosing this product from TFA Dostmann.

1. Before you use this product

- Please make sure you read the instruction manual carefully.
- Following and respecting the instructions in your manual will prevent damage to your instrument and loss of your statutory rights arising from defects due to incorrect use.
- We shall not be liable for any damage occurring as a result of non-following of these instructions. Likewise, we take no responsibility for any incorrect readings or for any consequences resulting from them.
- Please pay particular attention to the safety notices!
- Please keep this instruction manual safe for future reference.

2. Delivery contents

- Bathroom scale
- Batteries 2 x CR2032
- Instruction manual

3. Field of operation and all the benefits of your new instrument at a glance

- Automatic power-on function – no tapping required
- Automatic power-off function
- Modern flat design
- Tempered safety-glass surface

4. For your safety

- This product is exclusively intended for the field of application described above. It should only be used as described within these instructions.
- Before starting a treatment or diet discuss with your doctor first.

! Caution!
Risk of injury:

- Do not stand on the scale with wet feet and especially not when the surface of the scale is wet. There is a dangerous risk of slipping.
- Do not step on the edge of the scale: Danger of tilting!
- Keep this device and the batteries out of reach of children.
- Small parts can be swallowed by children (under three years old).
- Batteries contain harmful acids and may be hazardous if swallowed. If a battery is swallowed, this can lead to serious internal burns and death within two hours. If you suspect a battery could have been swallowed or otherwise caught in the body, seek medical help immediately.
- Batteries must not be thrown into a fire, short-circuited, taken apart or recharged. **Risk of explosion!**
- Low batteries should be changed as soon as possible to prevent damage caused by leaking.
- Never use a combination of old and new batteries together, nor batteries of different types.
- Avoid contact with skin, eyes and mucous membranes when handling leaking batteries. In case of contact, immediately rinse the affected areas with water and consult a doctor.

! Important information on product safety!

- Do not expose the device to extreme temperatures, vibrations or shocks.
- Do not place any objects on the unit.
- Protect it from humidity.

5. Set up and operation

- Open the battery compartment at the back and remove the insulation strip. Close the battery compartment again.
- All LCD segments will light up. "0.0 kg" (default: kg) appears on the display.
- Your scale is now ready to use.
- Press the UNIT button at the back and select your desired weight unit.
- Place the scale on a flat solid surface.
- In order to switch on the scale, just stand on the platform.
- After a few seconds your measured weight will be displayed in kg (default) or in the weighting unit you have set.
- The measuring process is finished as soon as your weight flashes on the display.
- The weight remains on the display after leaving the scale.
- The scale will turn off automatically after 10 seconds.

5.1 Setting of the weight unit

- At the back of the scale you will find the UNIT button to convert the weight unit. You can select between "kg" for Kilogram, "lb" for British and American pound (1 pound = 453,59237 Gram) and "st" for British Stone (1 st = 6,35 kg).

5.2 Use of the scale

- If possible, always weigh yourself at the same time of day (ideally in the morning), after going to the toilet, on an empty stomach and without clothes to achieve comparable results.

Instruction manual
Art.-No. 50.1015.01

6. Care and maintenance

- Clean your instrument with a soft damp cloth. Do not use solvents or scouring agents.
- Remove the batteries if you do not use it for a long period of time.
- Keep your scale in a dry place.

7. Troubleshooting

Problems	Solution
No indication	→ Stand on the platform or press the UNIT button → Ensure that the batteries' polarity is correct (+ pole above) → Change the batteries
Indication "Err"	→ Overload (max. 200 kg)
Indication "Lo"	→ Change the batteries
Incorrect display	→ Change the batteries

If your device fails to work despite these measures contact the supplier from whom you purchased it.

8. Waste disposal

This product and its packaging have been manufactured using high-grade materials and components which can be recycled and reused. This reduces waste and protects the environment.

Dispose of the packaging in an environmentally friendly manner using the collection systems that have been set up.

Disposal of the electrical device

Remove non-permanently installed batteries and rechargeable batteries from the device and dispose of them separately.

This product is labelled in accordance with the EU Waste Electrical and Electronic Equipment Directive (WEEE). This product must not be disposed of in ordinary household waste. As a consumer, you are required to take end-of-life devices to a designated collection point for the disposal of electrical and electronic equipment, in order to ensure environmentally-compatible disposal. The return service is free of charge. Observe the current regulations in place!

Disposal of the batteries

Batteries and rechargeable batteries must never be disposed of with household waste. They contain pollutants such as heavy metals, which can be harmful to the environment and human health if disposed of improperly, and valuable raw materials such as iron, zinc, manganese or nickel that can be recovered from waste. As a consumer, you are legally obliged to hand in used batteries and rechargeable batteries for environmentally friendly disposal at retailers or appropriate collection points in accordance with national or local regulations. The return service is free of charge. You can obtain addresses of suitable collection points from your city council or local authority.

The names for the heavy metals contained are:

Cd=cadmium, Hg=mercury, Pb=lead.

Reduce the generation of waste from batteries by using batteries with a longer lifespan or suitable rechargeable batteries. Avoid littering the environment and do not leave batteries or battery-containing electrical and electronic devices lying around carelessly. The separate collection and recycling of batteries and rechargeable batteries make an important contribution to relieving the impact on the environment and avoiding health risks.

! WARNING!
Damage to the environment and health through incorrect disposal of the batteries!

! WARNING!
Batteries containing lithium can explode

- Batteries and rechargeable batteries containing lithium (Li=lithium) present a high risk of fire and explosion due to heat or mechanical damage with potentially serious consequences for people and the environment. Pay particular attention to correct disposal.

9. Specifications

Measuring range	5 – 200 kg / 440 lb / 31 st
Resolution	0,1 kg / 0,2 lb
Batteries	2 x button cell batteries CR2032 (included)
Brand name	Tianqiu
Manufacturer	TianQiu Industrial Park, TianQiu Industrial Park, Xinji Industrial Zone, Machong Town, Dongguan, GuangDon, China
Dimensions	320 x 320 x 22 mm
Weight	1690 g (instrument only)

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No part of this manual may be reproduced without written consent of TFA Dostmann. The technical data are correct at the time of going to print and may change without prior notice. The latest technical data and information about this product can be found in our homepage by simply entering the product number in the search box.

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TFA Dostmann GmbH & Co.KG